



THE LAB EQUIPMENT

Email : [sales@thelabequipment.com](mailto:sales@thelabequipment.com)

Phone: +91-8759245890

**Product Name :**

Venturi Face Mask Adult Single Use

**Product Code :**

SLE/RAE/16



**Website:** <https://www.thelabequipment.com>, **Email:** [sales@thelabequipment.com](mailto:sales@thelabequipment.com)

plot. no.296. Sector 5. Vaishali. Ghaziabad. Uttar Pradesh 201010



THE LAB EQUIPMENT

Email : [sales@thelabequipment.com](mailto:sales@thelabequipment.com)

Phone: +91-8759245890

## Description :

### DESCRIPTION

Face mask to provide oxygen at flow rates between 1 – 10 LPM. Allows for precise oxygen delivery by utilizing different sized port to change FiO<sub>2</sub> (24 – 50%). Also known as air-entrainment masks able to provide total inspiratory flow at a specified FiO<sub>2</sub>.

### SPECIFICATIONS

**Material:** mask and tubing made from medical grade PVC

The mask has an adjustable nose clip and an adjustable strap to ensure the mask can be fitted securely.

**Length of the oxygen tube:** 1.5 to 2.1m

It delivers oxygen, with a specific concentration from 24% to 50% minimum "

The kit of the mask includes the tubing, humidity cup and multiple jets, which are color-coded and indicating the percentage of oxygen.

Single use, non-sterile item.

Adult size

FiO<sub>2</sub> ranges 24% (4 LPM), 28 % (4 LPM), 31% (6 LPM), 35% (8 LPM), 40% (8 LPM), 50% (10 LPM)

Anti-crush tubing is terminated with a "standard" connector

### ENVIRONMENTAL CONDITIONS

**Storage conditions:** 0 to 50 deg. C, 10 to 90% relative humidity, non-condensing

**Operating conditions:** 5 to 40 deg. C, 10 to 90% relative humidity, non-condensing

Contact Science Lab for your Educational School Science Lab Equipments. We are best technical educational equipments manufacturers, technical educational equipments manufacturer, technical educational equipments supplier, technical educational instruments, technical educational instruments exporter, technical educational instruments manufacturer, vocational training lab instruments exporter.